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SUMMARY

Silver Knight™ is an anesthesia breathing system tubing with an anti-microbial filter that uses silver ions to permit its use for a period of seven days. We (UETMIS) were asked to evaluate if this system should replace the single-use circuits currently used in our hospital.

No evidence-based data from controlled randomised or prospective studies regarding the effectiveness and safety of this system when used over seven days were found.

The analysis of data from national and international government agencies as well as those published in specialized journals using methodological and ethical validity criteria leans towards the principle of caution. We conclude that a single-use filter placed on the external breathing circuit between the patient and the anesthesia machine is required to reduce the risk of microbial contamination. Given our institution's goal of carbon neutrality in 2040 we suggest proceeding with a well controlled study at our center to evaluate the safety of Silver Knight (or similar) technology.

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AAGBI	Association of Anaesthetists of Great Britain and Ireland
AETMIS	Agence d'évaluation des technologies et des modes d'intervention en santé
ANZCA	Australian and New Zealand College of Anaesthetists
ASA	American Society of Anesthesiologists
BIS	British and International Standards
CDC	Centers for disease control and prevention
DGAI	German Society for Anaesthesiology and Intensive Care
ECR	Essai clinique randomisé
ETMISSS	Évaluation des technologies et des modes d'intervention en santé et services sociaux
FE	Filtre électrostatique
FHP	Filtre hydrophobe plissé
INSPQ	Institut national de santé publique du Québec
MHRA	Medicines and Healthcare Products Regulatory Agency
NaCl	Chlorure de sodium
PCI	Prévention et contrôle des infections
SARM	<i>Staphylococcus aureus</i> résistant à la mêtécilline

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¹ Assessing breathing-system filters. AR Wilkes - Medical device technology, 2004 - europepmc.org.

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6 RÉFÉRENCES

1. Centers for Disease Control and Prevention. Guidelines for preventing health-care-associated pneumonia, recommendations of CDC and the Healthcare Infection Control Practices Advisory Committee. MMWR 2004;53(No. RR-3):1-36
2. Klompas, M., et al., Strategies to prevent ventilator-associated pneumonia in acute care hospitals: 2014 update. Infect Control Hosp Epidemiol, 2014. 35 Suppl 2: p. S133-54.

⁴ <https://www.intersurgical.com/products/anaesthesia/modular-flexitube-silver-knight-breathing-systems>.

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Date de la recherche : 01-01-2000 au 01-10-2023

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Pubmed

#1 "Filtration/instrumentation"[Mesh] OR filter*

#2 "Anesthesiology"[Mesh] OR "Anesthesia"[Mesh] OR anesthe* OR anaesthe*

#3 "Cross Infection"[Mesh] OR "Equipment Contamination"[Mesh] OR infection* OR contamination* OR "Respiratory

Tract Infections"[Mesh]

#4 #1 AND #2 AND #3

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Cochrane Library

#1 MeSH descriptor: [Filtration] explode all trees and with qualifier(s): [Instrumentation - IS]

#2 filter*

#3 #1 or #2

#4 MeSH descriptor: [Anesthesia] explode all trees

#5 MeSH descriptor: [Anesthesiology] explode all trees

#6 anesthe* or anaesthe*

#7 #4 or #5 or #6

#8 MeSH descriptor: [Cross Infection] explode all trees

#9 MeSH descriptor: [Equipment Contamination] explode all trees

#10 MeSH descriptor: [Respiratory Tract Infections] explode all trees

#11 infection* or contamination*

#12 #8 or #9 or #10 or #11

#13 #3 and #7 and #12

Limites : Reviews, Technology Assessments

Embase

#1 'breathing circuit bacterial filter'/exp OR filter*

#2 'anesthesia'/exp OR anesthe* OR anaesthe*

#3 'hospital infection'/exp OR 'contamination'/exp OR 'respiratory tract infection'/exp OR infection* OR contamination*

#4 #1 AND #2 AND #3

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- PROSPERO Center for Reviews and Dissemination <http://www.crd.york.ac.uk/prospero/>
- Cochrane The Cochrane Library www.thecochranelibrary.com
- U.S. National Institute for Health Research <http://www.Clinicaltrials.gov>
- Current Controlled Trials Ltd. <http://www.controlled-trials.com>

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To : Dr. Lyndia Dervis, Aved Haroutounian

From : Freed-Lynn Noël, Rosalba Guarnieri, Lucia Di Cesare, Dr. Trong Tien Nguyen

Date : May 25, 2023

Subject : Re-use of Silver Knight anti-microbial breathing systems in the Operating Room at St. Mary's Hospital Center (Pilot project)

Infection Prevention & Control Service Summary

In the context of the operating room waste reduction projects led by Dr. Dervis, a risk evaluation for product re-use of the Silver Knight anti-microbial breathing system was requested from the Infection Prevention & Control Service. A thorough literature review of nosocomial risk associated with anesthesia devices was performed.

Infection Prevention & Control Service recommendations were provided.

Based on the manufacturer's claim that the Silver Knight breathing system could be used for up to 7 consecutive days, the device was implemented by the anesthesia and respiratory therapy service in February 2023 in operating room theaters 1, 3, 4, 7 & 8.

As part of this pilot project, a 4-week MRSA surveillance program was put in place to better understand the quality of disinfection. Following the 7-day use of the anti-microbial breathing system, and before disposing of the device, the OR was requested to obtain two cultures: a specimen culture of the exposed outer surface of the tubing and a culture of the inner surface of the tubing at the patient connection site.

During the 4-week surveillance, from April 9 to May 8, 2023 the laboratory service received 34 culture specimens from operating room theaters 1,3,4,7 & 8. All cultures results were negative for MRSA.

The surveillance showed no evidence of microbial growth, specifically MRSA, following the use of the 7- day Silver Knight anti-microbial breathing system in the operating room.

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PCI Recommendations

For continuous 7-day use of the silver Knight breathing system device in the operating room :

- Follow manufacturer instruction for use of the 7-day silver Knight breathing system.
- Follow Infection Prevention & Control Service recommendations sent December 9, 2022.
- A formal procedure must be written for the use of silver Knight breathing system before extending this pilot project to other CIUSSS sites.
- Notify PCI service of any changes in practice related to this device.
- PCI team to present the anesthesia pilot project to Infection Control Executives for review and approval.

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Intersurgical Silver Knight Breathing System for 7 Days Use (Intersurgical, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ)

Summary

Testing has been carried out to confirm the acceptability of using Intersurgical Silver Knight breathing systems for 7 days. Testing included bacterial reduction over 24 hours and 7 days, silver leaching over 24 hours and 7 days and testing according to ISO standards for anaesthetic and ventilator breathing systems. Microbial reduction was found to be effective after 7 days with no detectable leaching of silver from the tubes and the breathing systems were found to conform to ISO requirements after 7 days usage.

Intersurgical Silver Knight breathing systems are considered acceptable for 7 days use.

Method

7 days conditioning

The test for a single heated wire system (Ref: 2025000) was set up as in figure 1. Ventilator breathing settings were 500 ml TV, 20 breaths per minute and 1:2 I:E ratio.

Figure 1 - Heated Wire Experimental setup



Figure 2 - Uniflow Experimental setup

A non-heated wire system (Ref: 2000000) was set up equivalently but with no humidifier in place and an HMEF at the patient connection. The HMEF was replaced daily in accordance with accepted clinical practice.

The test for a UniFlow system (Ref: 2900000) was set up as in figure 2 with Patient breathing settings of 600 ml TV, 15 breaths per minute, 1:2 I:E ratio and 2L/min fresh gas flow.

Bacterial Reduction

The Nosocomial Infection National Surveillance Service (NINSS) report on the Surveillance of Hospital Acquired Bacteraemia indicates that there are a number of microorganisms that contribute to the majority of Hospital Acquired Infections. The following range of organisms were evaluated to establish the extent of effectiveness of the anti-microbial tubing :

- Methicillin Resistant Staphylococcus Aureus (KB-1005)
- Staphylococcus epidermis (NBRC-12689)
- Pseudomonas aeruginosa (IFO-13277)
- Klebsiella pneumoniae (NBRC-13889)
- Acinetobacter calcoaceticus (IFO-13006)
- Escherichia coli

■ **Preparation of bacterial inoculum**

Nutrient broth was diluted 500 times and adjusted to pH7. This 1/500 nutrient broth media was transplanted with the bacterial strain and used as the inoculum.

■ **Preparation of test samples**

Samples of 50 mm tube were cut from the Silver Knight product. All surfaces of the test pieces were cleaned with gauze containing ethanol and allowed to dry.

■ **Procedure**

Each test piece was placed in sterilised petri dishes and 0.5ml of the inoculum (containing $>1.0 \times 10^5$ colonies) was dropped onto the surface of the test piece. The sample was then covered with the lid of the petri dish and incubated at $35 \pm 1^\circ\text{C}$ and RH 90 % or higher for the test period. A test control and a control breathing system containing no Silver Knight anti-microbial additive were run in duplicate to ensure strain viability and test consistency.

■ **Measurement of viable cells**

At the end of the test time the inoculum was washed off the samples with 10 ml of diluent. Following serial dilution on to Tryptone Soya Agar and incubation at $35 \pm 1^\circ\text{C}$ for 40-48 hours the remaining viable cells were counted and the log reduction in growth established.

■ **E.Coli method**

Silver Knight respiratory tubing underwent 7 day simulated use conditioning. Samples were then taken and placed in individual petri dishes. Antibacterial activity was determined using the method outlined in ISO 22196:2007.

Silver Leaching

Cobalt (internal standard) was added to an acidified (HNO_3) known weight of sample. Silver was determined in the prepared solutions by Inductively Coupled Plasma Optical Emission Spectroscopy (ICP) against certified standards.

Results

Resistance to flow of the Uniflow breathing system was measured according to the method described in annex A of ISO 8835-2:2007. Resistance to flow for the heated and non-heated wire breathing systems was measured according to the method described in ISO 5367:2000.

- **Leakage** : Leakage testing was carried out according to ISO 5367:2000.
- **Tapers** : Taper testing was carried out according to ISO 53561:2004

Table 1 - Silver leaching from Silver Knight breathing system over 24 hours and 168 hours

Reference	Silver mg/kg
Limit of Detection	0.01
CONTROL	<0.01
24 hours 1	<0.01
24 hours 2	<0.01
7 days 1	<0.01
7 days 2	<0.01

Table 2 - Resistance to flow measurements for the Uniflow system according to ISO 8835-2:2007 after 7 days usage

SAMPLE	MAX PRESSURE (MBAR)	MIN PRESSURE (MBAR)
2900000	4.2	2.1

Table 3 - Resistance to flow for the 2000000 and 2025000 systems at 60 L/min after 7 days usage

SAMPLE	RESISTANCE TO FLOW (MBAR)	
	Inspiratory	Expiratory
2000000	0.6	0.6
2025000	1.0	1.0

Table 4 - Leakage for the 2000000, 2025000 and 2900000 systems at 60 mbar after 7 days usage

SAMPLE	LEAKAGE (ML/MIN)		
	2000000	2025000	2900000
1	<5ml/min	<5ml/min	<5ml/min
2	<5ml/min	<5ml/min	<5ml/min

Table 5 – Taper sizes for all connections on the Uniflow system after 7 days usage

CONNECTION	TAPERS
Patient End 22M	Small
Patient End 15F	Medium
Machine End 22F	Medium-Small
Machine End 22M	Small

Table 6 – Taper sizes for connections on 2000000 and 2025000 systems after 7 days usage

CONNECTION	TAPERS	
	2000000	2025000
Patient End 22M	Medium	Medium
Patient End 22F	Small	Small
Patient End 22F	Small	Small

Discussion

Microbial Reduction

Graph 1 shows that there is a significant reduction in the number of viable cells on the Flextube containing Silver Knight additive compared to the control Flextube samples with no Silver Knight additive. This reduction is seen to be consistent after both 24 hours and 7 days.

Silver Leaching

Table 1 shows that there was no detectable silver additive found to leach from the tubing after 24 hours or 7 days.

Resistance to flow, Leakage and Tapers

Comparing the results from tables 2-6 to the requirements from the relevant standards, the results show that all the breathing systems conform to the requirements of each standard after 7 days usage.

For Resistance to Flow, Compliance, Leakage and Tapers, the testing has been carried out on standard breathing systems without the Silver Knight additive in the Flextube. The only difference between standard breathing systems and Silver Knight systems is the anti-microbial additive and as such testing on the standard systems can be considered equivalent to testing on Silver Knight breathing systems.

Conclusion

Intersurgical Silver Knight breathing systems are considered acceptable for 7 days use.

References

- NINSS reports on surgical site infection and hospital acquired bacteraemia. Commun Dis Rep CDR Wkly. 2000 Jun 16;10(24):213, 216
- ISO 22196:2007
- ISO 8835-2:2007
- ISO 5367:2004
- ISO5356-1-2004
- Silver Knight Leaching Results
- Critical care systems for 7 days use
- JANN_44848_Anaes_UniFlow 7 day validation_issue1

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